## Airmaxx CAM Walker Boot (High Top)

Suggested HCPCS: L4360 / L4361





## **Indications**

Used for soft tissue injuries, grade 2 and 3 sprains stable fractures and post-operative stabilization. Helps to reduce pain, swelling and edema of foot and ankle.



- Inflatable bladder incorporated into liner for increased immobilization of the lower leg, ankle and foot.
- Specially designed reinforced adjustable toe guard for added protection.
- Cushioned sole for greater energy absorption.
- State of the art design provides greater protection without sacrificing comfort.
- Fits either left or right fool.

## **Sizes**

Item#	Size	Male Shoe Size	Female Shoe Size
AB0004-141-00	X-Small	Up to 4	Up to 5 <sup>1/2</sup>
AB0004-141-01	Small	41/2 - 7	6-8
AB0004-141-02	Medium	71/2 - 10	81/2 - 111/2
AB0004-141-03	Large	101/2 - 121/2	111/2 - 131/2
AB0004-141-04	X-Large	121/2+	131/2+

P: 800.216.1479 | F: 888.335.5394